

Quantification of anti-mTNF-alpha ADCP activity using *iLite*[®] ADCP Effector Assay Ready Cells

For research and professional use only. Not for use in diagnostic procedures.

This application note contains a suggested protocol and performance data. Each individual laboratory must set up their own method and perform relevant validations.

Background

The immune system uses various mechanisms to kill specific pathogens, infected cells, and cancer cells. Therapeutic antibodies act by binding to a cell surface receptor by the Fab domain resulting in induction/blocking of signaling events but also the Fc part of the antibody has an important biological role. Fc engineering strategies to increase the efficacy of anti-cancer antibodies are ongoing (1). The Fc-part of the antibody are involved in inducing antibody-dependent cellular cytotoxicity (ADCC), antibody-dependent cellular phagocytosis (ADCP), complement-dependent cytotoxicity (CDC), and programmed cell death (1).

ADCP is a mechanism whereby pathogenic cells are phagocytosed by various immune effector cells such as monocytes, macrophages, dendritic cells, and neutrophils (2). The antibody binds to a specific antigen on the pathogen/cancer cell, while the Fc part of the antibody binds and activates the Fc receptor on the effector cell resulting in activation of a pathway that triggers phagocytosis. The most dominant receptor in induction of ADCP is the low affinity Fc receptor FcγRIIa (CD32a).

TNF-alpha is expressed as a membrane bound trimeric type II transmembrane protein which is cleaved by the metalloprotease TNF-converting enzyme (ADAM17) (3). Due to diverse stimuli the soluble form of TNF-alpha is generated (4). The first “biologic” on the market was a monoclonal antibody, Infliximab, against TNF-alpha and it was used in rheumatoid arthritis (5). The mode of action of anti-TNF antibodies are generally used to block TNF-alpha activities but the therapeutic efficacy of the anti-TNF antibodies have been shown to depend on the Fc-Fcγ receptor interaction of the antibodies (6).

The use of the *iLite*[®] ADCP Assay Ready Cells together with target cells of interest offers a technology platform to study the ADCP mechanism of a therapeutic anti-TNF alpha antibody in a reliable and reproducible manner.

Principle of the assay

The *iLite*[®] ADCP Assay Ready Cells are engineered cells that enable ADCP to be examined through the specific expression of Firefly luciferase. When the antibodies of interest bind to the antigens on the surface of the target cell, the target-bound antibodies will be presented to FcγRIIa on the effector cells. When the Fc-portion of the target-bound antibodies binds to the receptor, multiple cross-linking of the two cell types occurs. This will initiate a signaling cascade through the NFAT response element which triggers the expression of Firefly luciferase (FL) in the effector cells. In this application note, we describe the use of an effector cell line (*iLite*[®] ADCP Effector Assay Ready Cells) that over-express FcγRIIa and contain the FL reporter gene that responds to the principal transcription factors that mediate signaling from the FcγRIIa receptor, together with a positive target cell line which expresses membrane bound TNF-alpha (*iLite*[®] mTNF-alpha (+) Target Assay Ready Cells). *iLite*[®] ADCP Effector Assay Ready Cells also contain the Renilla Luciferase (RL) reporter gene, under the control of a constitutive promoter, that allows drug-induced FL activity to be normalized with respect to the constitutive expression of RL.

This render assay results independent of variations in cell number, serum matrix effects, or lysis of the effector cells by the target cells. In addition, we also describe the use of a negative control in the form of a target cell line without membrane bound TNF-alpha expression (*iLite*[®] mTNF-alpha (-) Target Assay Ready Cells). The Firefly luciferase signal can be measured in a luminometer following addition and incubation of luciferase substrate. The Firefly luciferase signal is proportional to the functional activity of anti-TNF antibodies such as Infliximab in the sample (Fig.1).

Specimen collection

Quantification of anti-mTNF-alpha ADCP activity using *iLite*[®] ADCP Effector Assay Ready Cells together with *iLite*[®] mTNF-alpha Target Assay Ready Cells can be performed in test samples including human serum.

Material and equipment needed

Material and equipment	Suggested supplier	Reference
<i>iLite</i> [®] ADCP Effector Assay Ready Cells	Svar Life Science	BM5004
<i>iLite</i> [®] mTNF-alpha (+) Target Assay Ready Cells	Svar Life Science	BM5013
<i>iLite</i> [®] mTNF-alpha (-) Target Assay Ready Cells	Svar Life Science	BM5014
Diluent (RPMI 1640 + 9% heat inactivated FBS + 1% Penicillin Streptomycin)	Gibco	61870 (RPMI) 26140-079 (FBS) 15140-122 (Penicillin-Streptomycin)
Infliximab or analogues	Janssen Biotech	NA
Firefly/Renilla luciferase substrate	Promega	E2940, Dual-Glo [®] Luciferase Assay System
Plate; White walled micro well plate suitable for luminescence	PerkinElmer	6005680
Microplate Luminometer with appropriate reading software – no filter on luminometer	Contact Svar Life Science for list of recommended suppliers	NA
Incubator, 37°C with 5% CO ₂	NA	NA
Water bath, 37°C	NA	NA
Single-channel and multi-channel pipettes with polypropylene disposable tips	NA	NA
Polypropylene tubes or plate for dilution	NA	NA
Single-use polypropylene reservoir	NA	NA
Plate shaker	NA	NA
Timer	NA	NA

Preparation of calibrators (example given with infliximab)

The ADCP effect of the Infliximab antibody from Janssen Biotech has successfully been measured in combination with a mix of ADCP Effector Assay Ready Cells and mTNF-alpha (+) Target Assay Ready Cells. As a negative control, a combination of ADCP Effector Assay Ready Cells and mTNF-alpha (-) Target Assay Ready Cells was used. In the present assay an Effector:Target ratio of 4:3 has been used. The optimal ratio is dependent on the antibody and target cells used and should be determined each

time a new assay is set up. The table below shows recommended dilutions of Infliximab when making an 8-point calibration curve.




Figure 1. Example of Infliximab calibration curve using Firefly Luciferase substrate Dual-Glo® Luciferase Reagent. Values are shown as mean of triplicate \pm SD and values on x-axis are given as **final concentration** in the wells before addition of Dual-Glo® Luciferase Reagent.

Calibrator	Infliximab	
	Calibrator	solution conc. (ng/mL)
1		1 000
2		500
3		200
4		100
5		50
6		25
7		6.0
8		0

Table 1. Suggested calibrator **solution concentrations** for Infliximab.

Protocol

Assay preparation and incubation

1. Design a plate layout.
2. Dilute calibrators, controls, and samples to fall within the expected **in assay values (= final concentration)** of 0-500 ng/mL.
3. Add 40 μ L calibrators, controls, and samples in duplicate to assigned wells.
4. Thaw the vial of ADCP Effector Assay Ready Cells and the vials of mTNF-alpha (+) Target Assay Ready Cells and mTNF-alpha (-) Target Assay Ready Cells in a 37°C water bath with gentle agitation.
5. Mix the cell suspensions very carefully **at least 10 times with a pipette** in order to ensure a homogeneous distribution of cells.
6. Dilute 200 μ L of the ADCP Effector Assay Ready Cells and 200 μ L the mTNF-alpha (+) Target Assay Ready Cells with 3.44 mL Diluent. The total volume of the diluted ADCP Effector /Target mTNF-alpha (+) Assay Ready Cells mixture is 3.84 mL.
7. In a separate tube, dilute 50 μ L of the ADCC Effector (V) Assay Ready Cells with 50 μ L of the mTNF-alpha (-) Target Assay Ready Cells with 860 μ L Diluent. The total volume of the diluted ADCC Effector (V)/Target mTNF-alpha (-) Assay Ready Cells mixture is 960 μ L.
8. Add 40 μ L of the diluted cells to each well to be tested.
9. Place the lid on the plate and mix on a plate shaker at **minimum of 750 rpm** for 10 sec. Alternatively, mix the cell suspensions very carefully in the wells by pipette. Insufficient mixing can cause reduced assay sensitivity.
10. Incubate for 4 hours at 37°C with 5% CO₂.

Adding substrate solutions

11. Equilibrate the plate and the substrate solutions to room temperature.
12. Prepare the **Firefly luciferase substrate** in accordance with the supplier's instructions and add 80 μ L per well. Mix and protect the plate from light. After 10 minutes incubation at room temperature read in a luminometer.
13. If appropriate, prepare the **Renilla luciferase substrate** in accordance with the supplier's instructions and add 80 μ L per well. Mix and protect the plate from light. After 10 minutes incubation at room temperature read in a luminometer.

Normalization

The reporter gene used for normalization, Renilla luciferase, is under the control of a tyrosine kinase promoter, and is thus constitutively expressed. Unspecific effects such as serum matrix effects or differences in cell number can be obviated by relating the specific Firefly signal with the Renilla normalization signal through simple division.




Figure 2.

A: Measurement of the specific Firefly (FL) signal. Cells were stimulated with increasing concentrations of Infliximab.

B: Measurement of the Renilla (RL) signal for normalization, from the identical assay as in A.

C: Dose response curve after normalization of the specific Firefly signal with the according Renilla signal.

Precautions


- This application note is intended for professional laboratory research use only. The data and results originating from following the Application Note should not be used either in diagnostic procedures or in human therapeutic applications.
- Use and handle the material and instruments referenced according to the suppliers'/manufacturers' instructions or product specifications accompanying the individual material and instruments.
- Dispose of all sample specimens, infected or potentially infected material in accordance with good microbiological practice. All such materials should be handled and disposed as though potentially infectious.
- Residues of chemicals and preparations are generally considered as biohazardous waste and should be inactivated prior to disposal by autoclaving or using bleach. All such materials should be disposed of in accordance with established safety procedures.

Proprietary Information

In accepting delivery of *iLite*[®] Assay Ready Cells the recipient agrees not to sub-culture these cells, attempt to sub-culture them or to give them to a third-party recipient, and only to use them directly in assays. *iLite*[®] cell-based products are covered by patents which are the property of Svar Life Science AB and any attempt to reproduce the delivered *iLite*[®] Assay Ready Cells is an infringement of these patents.

QUICK GUIDE

Quantification of anti-mTNF-alpha ADCP activity using *iLite*[®] ADCP Effector Assay Ready Cells



Troubleshooting and FAQ

Please consult the Svar Life Science website www.svarlifescience.com

References

1. Liu R, Oldham RJ, Teal E, Beers SA, and Cragg MS. *Fc-Engineering for Modulated Effector Functions—Improving Antibodies for Cancer Treatment*, *Antibodies*. 9(4):64 (2020)
2. Graziano RF and Engelhardt JJ. *Role of FcγRs in Antibody-Based Cancer Therapy*, *Curr Top Microbiol Immunol*. 423:13-34 (2019).
3. Aggarwal, B. B., Gupta, S. C. & Kim, J. H. *Historical perspectives on tumor necrosis factor and its superfamily: 25 years later, a golden journey*. *Blood* 119, 651–665 (2012).
4. Black, R. A. et al. *A metalloproteinase disintegrin that releases tumour-necrosis factor-α from cells*. *Nature* 385, 729–733 (1997)

5. Monaco C, Nanchahal J, Taylor P and Feldmann M. Anti-TNF therapy: past, presence and future. *Int. Immunol.*27(1):55-62 (2015)
6. McRae BL, Levin AD, Wildenberg ME, et al. *Fc Receptor-mediated effector function contributes to the therapeutic response of anti-TNF monoclonal antibodies in a mouse model of inflammatory bowel disease.* *J Crohns Colitis* 2016;10:69–76.