

Quantification of functional IFN alpha using *iLite*[®] IFNI FAST Assay Ready Cells

For research and professional use only. Not for use in diagnostic procedures.

*This application note contains a suggested protocol and performance data.
Each individual laboratory must set up their own method and perform relevant validations.*

Background

Interferons (IFNs) are cytokines with antiviral, antitumor and immunoregulatory functions, released by both immune and non-immune cells as a first line of host innate defense, particularly viral infections. IFNs are classified into three groups and generally initiate signaling via the JAK-STAT (Janus kinase/signal transducer and activator of transcription) pathway. Type I IFNs includes at least 16 different subtypes, most well defined are IFN alpha (IFN α) and IFN beta (IFN β) which binds to the heterodimeric IFN alpha/beta receptor (IFNAR) thereby inducing intracellular signaling (1,2). Type I IFNs are used in several clinical settings, due to their antiviral and immunomodulatory properties (3). These properties have generated interest in their clinical use to enhance antigen presentation, control viral infections and promote antitumor responses (3–5). Recently, a robust type I interferon response has been observed in patients with severe COVID-19, which could contribute to the exacerbated hyperinflammation observed in the progression to severe COVID-19 (6). Specifically, the administration of IFN alpha2 as monotherapy or in combination with other treatments resulted in tumor regression and/or prolonged survival in a diversity of other hematological and disseminated solid malignancies including myeloma, lymphomas, myeloproliferative syndromes, melanoma, renal cell and bladder carcinomas and Kaposi sarcoma (5).

Principle of the assay

The *iLite*[®] IFNI FAST Assay Ready Cells are engineered cells, optimized to express Firefly luciferase (FL) under the control of a type I IFN responsive promoter. *iLite*[®] IFNI FAST Assay Ready Cells also contain the Renilla Luciferase (RL) reporter gene, under the control of a constitutive promoter. The constitutive expression of RL allows normalization of type I IFN induced FL activity, and renders assay results independent of variations in cell number or serum matrix effects. The luciferase signal can be measured in a luminometer following addition and incubation of luciferase substrate. The Firefly luciferase signal is proportional to the concentration of functional type I IFN in a sample, in this application note exemplified by IFN alpha 2b (Fig.1).

Specimen collection

The *iLite*[®] IFNI FAST Assay Ready Cells can be used for measuring concentration of IFN alpha in test samples including human serum.

Material and equipment needed

| Material and equipment | Suggested supplier | Reference |
|---|---|--|
| <i>iLite</i> [®] IFNI FAST Assay Ready Cells | Svar Life Science | BM4049 |
| Diluent (DMEM containing 9% heat inactivated FBS + 1% Penicillin-Streptomycin) | Gibco | 31966-021 (DMEM) 26140-079 (FBS) 15140-122 (Penicillin-Streptomycin) |
| Human IFN alpha 2b (used in example protocol) | Immunotools | 11343516 |
| Firefly/Renilla luciferase substrate | Promega | E2920, Dual-Glo Luciferase Assay System |
| Plate; White walled micro well plate suitable for luminescence | PerkinElmer | 6005680 |
| Microplate Luminometer with appropriate reading software – no filter on luminometer | Contact Svar Life Science for list of recommended suppliers | NA |
| Incubator, 37 °C with 5% CO ₂ | NA | NA |
| Water bath, 37 °C | NA | NA |
| Single-channel and multi-channel pipettes with polypropylene disposable tips | NA | NA |
| Polypropylene tubes or plate for dilution | NA | NA |
| Single-use polypropylene reservoir | NA | NA |
| Plate shaker | NA | NA |
| Timer | NA | NA |

Preparation of type I IFN calibrator (example protocol given with rec. human IFN alpha 2b, Immunotools 11343516)

Recombinant human IFN alpha 2b from Immunotools has successfully been used to stimulate the *iLite*[®] IFNI FAST Assay Ready Cells.




Figure 1. Example of IFN alpha 2b calibration curve.

| Calibrator | IFN alpha 2b |
|------------|---|
| | Suggested calibrator solution conc. (IU/ml) |
| A | 20000 |
| B | 5000 |
| C | 1250 |
| D | 417 |
| E | 139 |
| F | 28 |
| G | 5.6 |
| H | 0 |

Table 1. Suggested calibrator solution concentrations for IFN alpha 2b.

Assay preparation and incubation

1. Design a plate layout. It is recommended to perform the test at least in duplicates.
2. Dilute calibrators, controls, and samples to fall within the expected **in assay values** of 0 - 10000 IU/mL.
3. Add 40 µL calibrators, controls, and samples in duplicate to assigned wells (final concentration will be half of solution concentration).
4. Thaw the vial of *iLite*[®] IFNI FAST Assay Ready Cells in a 37°C water bath with gentle agitation. The cell suspension is mixed very carefully ten times with pipette to ensure a homogeneous distribution of cells.
5. Dilute 250 µL cell suspension with 5.75 mL Diluent.
6. Add 40 µL diluted cells to each well.
7. Place the lid on the plate, mix and incubate for 5 hours at 37 °C with 5% CO₂.

Adding substrate solutions

8. Equilibrate the plate and the substrate solution to room temperature.
9. Prepare the **Firefly luciferase** substrate according to the manufacturer's instructions and add 80 µL per well. Mix and protect the plate from light. After 10 minutes incubation at room temperature read in a luminometer.
10. If appropriate, prepare the **Renilla luciferase** substrate according to the manufacturer's instructions and add 80 µL per well. Mix and protect the plate from light. After 10 minutes incubation at room temperature read in a luminometer.

Precautions


- This application note is intended for professional laboratory research use only. The data and results originating from following the Application Note should not be used either in diagnostic procedures or in human therapeutic applications.
- Use and handle the material and instruments referenced according to the supplier's/ manufacturer's instructions or product specifications accompanying the individual material and instruments.
- Dispose of all sample specimens, infected or potentially infected material in accordance with good microbiological practice. All such materials should be handled and disposed as though potentially infectious.
- Residues of chemicals and preparations are generally considered as biohazardous waste and should be inactivated prior to disposal by autoclaving or using bleach. All such materials should be disposed of in accordance with established safety procedures.

Proprietary Information

In accepting delivery of *iLite*[®] Assay Ready Cells the recipient agrees not to sub-culture these cells, attempt to sub-culture them or to give them to a third-party recipient, and only to use them directly in assays. *iLite*[®] cell-based products are covered by patents which are the property of Svar Life Science AB and any attempt to reproduce the delivered *iLite*[®] Assay Ready Cells is an infringement of these patents.

QUICK GUIDE

Quantification of functional IFN alpha using *iLite*[®] IFNI FAST Assay Ready Cells



Troubleshooting and FAQ

Please consult the Svar Life Science website www.svarlifescience.com

References

1. Ivashkiv LB, Donlin LT. Regulation of type I interferon responses. *Nat Rev Immunol.* 2013;14(1):36–49.
2. Padilla CML de, Niewold TB. The type I interferons: Basic concepts and clinical relevance in immune-mediated inflammatory diseases. *Gene.* 2016;576(1):14–21.
3. Park A, Iwasaki A. Type I and Type III Interferons – Induction, Signaling, Evasion, and Application to Combat COVID-19. *Cell Host Microbe.* 2020;27(6):870–8.
4. Propper DJ, Balkwill FR. Harnessing cytokines and chemokines for cancer therapy. *Nat Rev Clin Oncol.* 2022;19(4):237–53.
5. Borden EC. Interferons α and β in cancer: therapeutic opportunities from new insights. *Nat Rev Drug Discov.* 2019;18(3):219–34.
6. Lee JS, Shin EC. The type I interferon response in COVID-19: implications for treatment. *Nat Rev Immunol.* 2020;20(10):585–6.