

***iLite*[®] AAV Responsive Reporter Assay Ready Cells**

REF: BM6100

For research use only. Not for use in diagnostic procedures.

DESCRIPTION

iLite[®] AAV Responsive Reporter Assay Ready Cells are based on the human embryonic kidney 293 (HEK293) cell line and have been genetically engineered and optimized to be responsive to the transduction of AAV vectors produced by *iLite*[®] AAV Packaging Assay Ready Cells. The transduction leads to the activation of a transactivator regulating the Firefly Luciferase reporter gene construct. The cells also contain the Renilla Luciferase internal standardization gene under the control of a constitutive promoter that renders assay results independent of cell number and provides a means for correcting for cytotoxic effects that may be encountered with some biological samples.

CONTENT

>250 µL of Assay Ready Cells suspended in recovery medium from Gibco (Cat. No 12648-010).

RECEIPT AND STORAGE

Upon receipt confirm that adequate dry-ice is present, and the cells are frozen. Immediately transfer to -80°C storage. Cells should be stored at -80°C or at lower temperature and are stable as supplied until the expiry date shown. Cells should be diluted and plated immediately after thawing.

BACKGROUND

Currently, the AAV product development and approval processes have multiple challenges that require the development of new tools (1). The *iLite*[®] AAV Responsive Assay Ready Cells together with *iLite*[®] AAV Packaging Assay Ready Cells of any serotype represent a functional, robust, and agile technology to facilitate the determination and titration of neutralizing antibodies against AAV vectors.

APPLICATION

The *iLite*[®] AAV Responsive Reporter Assay Ready Cells together with *iLite*[®] AAV Packaging Assay Ready Cells of any serotype can be used as a transduction inhibition (TI) assay for the determination and titration of neutralizing antibodies against AAV vectors.

Application Notes for the following assays are available:

- Determination of NAbs against AAV in human serum using *iLite*[®] AAV Packaging and *iLite*[®] AAV Responsive Reporter Assay Ready Cells (LABEL-DOC-00580)
- Titration of NAbs against AAV in human serum using *iLite*[®] AAV Packaging and *iLite*[®] AAV Responsive Reporter Assay Ready Cells (LABEL-DOC-0599)

RELATED PRODUCTS

REFERENCES

SYMBOLS ON LABEL

PRECAUTIONS

PROPRIETARY INFORMATION

REF Product name

BM6002	<i>iLite</i> [®] AAV2 Packaging Assay Ready Cells
BM6005	<i>iLite</i> [®] AAV5 Packaging Assay Ready Cells
BM6006	<i>iLite</i> [®] AAV6 Packaging Assay Ready Cells
BM6008	<i>iLite</i> [®] AAV8 Packaging Assay Ready Cells
BM6009	<i>iLite</i> [®] AAV9 Packaging Assay Ready Cells

1. Gupta V, Lourenço SP, Hidalgo IJ. Development of Gene Therapy Vectors: Remaining Challenges. J Pharm Sci. 2020



Lot number



Temperature limitation



Catalogue number



Biological risk



Use by



Manufacturer

For research use only. This product is intended for professional laboratory research use only. The data and results originating from using the product should not be used either in diagnostic procedures or in human therapeutic applications.

iLite[®] AAV Responsive Reporter Assay Ready Cells are a stable transfected cell line of human origin classified as a Class 1 Genetically Modified Microorganism. They should be handled in accordance with EU directive (2009/41/EC) and disposed of in a licensed contained-use facility in accordance with these regulations. When used in accordance with the manufacturer's product specification, the requirements of EC Directive 2009/41/EC on the contained-use of genetically modified microorganisms are deemed to have been met.

Residues of chemicals and preparations generally considered as biohazardous waste should be inactivated prior to disposal by autoclaving or using bleach. All such materials should be disposed of in accordance with established safety procedures.

In accepting delivery of *iLite*[®] Assay Ready Cells the recipient agrees not to sub-culture these cells, attempt to sub-culture them or to give them to a third party, and only to use them directly in assays. *iLite*[®] cell-based products are covered by patents which is the property of Svar Life Science AB and any attempt to reproduce the delivered *iLite*[®] Assay Ready Cells is an infringement of these patents.