

CALPROSMART

DeepLink Integration

System requirements

- Host PRO app
- iOS 9 or higher/ Android 4.2 or higher
- Mobil phone application only



CalproSmart is an easy to use self-test for calprotectin measurement, and with a DeepLink access it can easily be integrated into Patient Reported Outcome (PRO) apps.

Combining the subjective PRO questionnaires with CalproSmart self-test results has many benefits, which include:


- Easier to collect more data points in studies and clinical trials
- Symptom scores are aligned with calprotectin measurements
- The treating physicians get a more wholesome overview of the patient followed up through e-health and distant management

DeepLink integration


- Customizable app-in-app solution
- Seamless DeepLink integration
- No registration, storage, or residual data in the CalproSmart system
- Detailed technical documentation for integration

CalproSmart solution


- Step-by-step guided app
- Quick assay – 20 min from start to results
- Wide range 77-1500 mg/kg
- Validated against laboratory assay




Before use, the host and the CalproSmart DeepLink app needs to be integrated by host app developers



The user will be required to have both apps installed on their phone



Registration, login, and data storage are done exclusively in the host app and its back-end system.



To perform the test, the user will click the designated button for this in the host app.



Results (and the user) is transferred back to the host app once the test is completed.